


# A-12, adapter

## DESCRIPTION

A-12, adapter for work with tubes which are 12 mm in external diameter.



## CAT. NUMBER

BS-050102-IK

A-12 adapter

## ACCESSORIES



**DEN-1**

BS-050102-AAF

Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units (0 –  $180 \times 10^7$  cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



**DEN-1B**

BS-050104-AAF

Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units (0 –  $180 \times 10^7$  cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)