

HSP-16/250, platform

DESCRIPTION

Platform with tight fit clamps for flasks.

SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	250-300 ml
Number of clamps	16



CAT. NUMBER

BS-010167-MK

Platform HSP-16/250

ACCESSORIES



ES-20/80

Software included, without platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C

Software included, without platform
Orbital Shaker-Incubator with cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)