

R-24GC, Rotor

DESCRIPTION

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)

SPECIFICATIONS

Rotor type	Swing-out
Card dimensions (l × w)	53 × 74 mm
Capacity	24 gel cards
Max. speed	1500 rpm
Max. RCF: LMC-3000	280 × g
Max. RCF: LMC-4200R	280 × g
Gel card Manufacturers	Grifols, DiaMed, etc.
Autoclavable	+



CAT. NUMBER

BS-010208-VK

Rotor R-24GC

ACCESSORIES



LMC-3000
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-12
Without rotor
Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-4200R
Without rotor
Laboratory Refrigerated
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors