

# Ri-12/15, Rotor

## DESCRIPTION

Swing-out rotor for conical bottom centrifuge tube.

## SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	17 × 120 mm
Capacity	12
Tube's volume	15 ml
Max. speed	4200 rpm
Max. RCF	3160 × g
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



## CAT. NUMBER

BS-010221-KK

Rotor Ri-12/15

## ACCESSORIES



**LMC-56**  
Without rotor  
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**RR-U**  
BS-010208-UK  
rack for rotors

Rack for rotors



**BN-17/120**  
BS-010208-TK  
bucket set for rotor R-12/15

Spare bucket set for rotor R-12/15. Plastic buckets (12 pcs).