

Ri-6, rotor

DESCRIPTION

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.

SPECIFICATIONS

Tube holders material	Aluminium
Rotor type	Swing-out
Tube dimensions (Ø × length)	29 × 115 mm
Capacity	6
Tube's volume	50 ml
Max. speed	4200 rpm
Max. RCF	3160 × g
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



CAT. NUMBER

BS-010221-HK

Rotor Ri-6

ACCESSORIES



LMC-56
 Without rotor
 Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BI-25-6
 BS-010221-VK
 adapter set

For rotors R-6, R-6P, Ri-6, Ri-6P optional insert set (6 pcs) - BI-25-6 for 25 ml conical tubes



RR-U
 BS-010208-UK
 rack for rotors

Rack for rotors