

SPM, double-sided adhesive mat

DESCRIPTION

Convenient alternative to traditional steel holders, an easy way to fix tubes, plates, flasks and other laboratory wares on platforms for cultivation, incubation and mixing.

SPM is compatible with PP-4 platform, that fits both on PSU-10i orbital shaker and in ES-20 Shaker-Incubator.

Made of polyurethane with adhesive, simple to clean and durable, able to withstand up to 1,000 times placement/removal or 12 months of use. Additional information about temperature, working volume and speed limitations available in user manual.

SPECIFICATIONS

| | |
|----------------------------|-------------------------------------------------------------|
| Colour | transparent |
| Duration of use | up to 1,000 times placement/removal or 12 months |
| Temperature range | +4C° ... +80C° |
| Speed range | 0–400 rpm |
| Overall dimensions (W×D×H) | 210 × 210 × 3 mm (double sided 1.5 mm PU with PET adhesive) |



CAT. NUMBER

BS-010111-BK

SPM, Double-sided adhesive mat

ACCESSORIES



ES-20
Without platform
Orbital Shaker-Incubator

The ES-20 is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program

Built-in ...

[read more](#)



PP-4
BS-010108-BK
platform

Flat platform with non-slip silicone mat for Petri dishes, culture flasks, agglutination cards.



PSU-10i
Without platform
Orbital Shaker

Shaker PSU-10i provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5)

...

[read more](#)